COMMUNIQUÉ DE PRESSE

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N° d'enregistrement: 5160407 VAT No: GB 838 5592 82

Illustration :

A Safety Guide for Using Laboratory Heating Blocks

Contributors: Dr I Smellie*, Dr N Langerman**, Dr Kerry Elgie***, Martyn Fordham***

Heating blocks provide a safe, convenient and productive alternative to heating mantles and oil baths for heating round-bottomed flasks, tubes and vials. Used in combination with a standard hotplate stirrer, heating blocks such as the Asynt DrySyn® range have proved themselves with their ability to outperform the heat-conducting properties of oil baths. They pose a far lower fire risk and their use makes the clean-up of glassware much easier as there is no residual oil contamination on the outside of the flasks. In addition to accelerating your chemical reactions, heating blocks ensure a safer, cleaner, healthier working environment.

However, solid aluminium heating blocks can reach temperatures of over 250 °C and therefore care needs to be taken not to touch the units, or inserts, at any time unless you are certain that they are cool. The nature of these heating blocks is such that the equipment will retain the heat for some time and larger heating blocks that have been used at high temperatures could stay hot enough to cause a burn even after over half an hour has elapsed. Ideally users should use a thermocouple probe in the heating block to monitor temperature during the cooling phase. In the absence of a thermocouple probe it is recommended that any block that has been heated should be left for a minimum of one hour.

In the case of Asynt DrySyn® Heating Blocks, you could also use an additional DrySyn® Safety Heat Sticker which shows when the heating block has cooled sufficiently to touch without causing a hazard.

This document has been written to provide laboratory staff with an **IN BEST PRACTICE GUIDE** for the safe use of heating block systems.

All DrySyn® bases are supplied with adjustable feet to offer compatibility with any hotplate stirrer.

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Asynt Ltd