ASVM

Lighting the Way to Sustainable Chemistry: The Power of Flow and Photochemistry

www.asynt.com



Photochemistry as a Greener Application

The power of flow and photochemistry.

- We'll explores how the synergy of continuous flow systems and light-driven reactions is transforming chemical synthesis—offering; enhanced safety, process intensification and reduced environmental impact
- Key themes include the benefits of precise reaction control, and improved energy efficiency in photochemical processes.
- The session will highlight recent advancements and practical applications, demonstrating how these technologies contribute to a more sustainable future in chemical manufacturing.







How do we Initiate a Chemical Reaction?

Chemical reactions are generally activated by heat or a chemical reagent

- One of the key advantages of modern laboratory equipment is its easier to use some of the lesser used methods:
 - Electrochemistry which uses electrons
 - Photochemistry which uses light
- An increasing trend in recent years to use both techniques in chemistry, sometimes referred to as reagentless (or traceless) chemistry.







in

Photochemistry – How It Began

Great masters of Organic Synthesis

The application of photochemistry for the synthesis of organic compounds is not a new concept:



Notably **Eaton** 1964 on his cubane synthesis.

- 2
- **Corey** 1964 applied photochemistry and his retrosynthetic approach to carry out the synthesis of caryophyllenes.
- 3
- Wender's synthesis of cedrene, in 1981.





Increased use of Photochemistry

Photochemistry Citations

Photochemistry research has seen an increase from researchers in academia and industry for several reasons.

- The development of commercially available photochemistry reactors.
- Novel synthetic routes reducing the number of steps.
- The development of novel photoredox catalysts.
- The awareness of photochemistry as a greener chemistry.



Estimated trend in publications citing "photochemistry" from 2000 - 2025



Recent Advances towards Sustainable Chemistry?

Photochemistry has gained popularity as a "greener" synthetic technique in recent years.

- The impacts in the synthetic photochemistry is already being seen in two main areas:
 - Photons can be used as a traceless reagent to activate otherwise benign starting materials
 - Molecules can be designed in fewer steps enabling savings in resources, time, energy and waste.

- While photochemistry has the potential in the advancement of sustainable processing there are still some hurdles in achieving this:
 - Light source technology.
 - Reactor design.
 - Process understanding.



12 Principles of Green Chemistry

- Photochemistry aligns closely with the 12 Principles of Green Chemistry, offering a sustainable approach to chemical synthesis by utilizing light as an energy source.
- Our focus is on the following factors:
 - The development of new light source technology. Better selectivity, improved efficiency, temperature control.
 - The advantages of modern photochemical reactor design.

Higher photon flux.







Photochemistry Light Sources

- Photochemistry in the past has utilised high intensity mercury lamps to facilitate photocatalyzed reactions
 - These produce a broad spectrum of wavelengths (200-600 nm).
 - They produce a lot of heat which is transferred to the reaction.
 - They almost always require external cooling to maintain temperature.
- Advances in LED technology have opened-up more targeted photochemistry
 - LEDs have a narrow emission band (normally in 20 nm range). This allows more selective chemistry.
 - LEDs have a small irradiation window which allows them to be directed toward the reactor.
 - LEDs also give off a lot less heat, so external cooling is not always required.





Photochemistry Light Sources



Effects of Mercury Vapour Pressure on the Spectral Power Distribution

¹ Lamptech.co.uk – Mercury Vapour





Photochemistry Light Sources



Emission spectrum of common light sources used for Photochemical Applications¹

Cambié, D., et al., Chem. Rev., 2016, 116 (17), 10276–10341

in



Photochemistry Reactors

The problem with light (photons) is they lose their energy very quickly.

This is described by the **Bouguer-Lambert-Beer law**, the further the reaction is from the light the less power it sees.

• When using traditional flasks, the reaction is only irradiated near the wall giving non-uniform reaction conditions

When performing photochemical reactions in flow reactors (with i.d. of < 1 mm), full transmission of light is allowed to the reaction solution.

• The result in a more homogeneous irradiation of the reaction giving shorter reaction times and less side-products caused by over-irradiation (often seen in equivalent batch conditions).



Attenuation of light with the distance of irradiation^{2,3}

Sambiagio, C.; Noël, T. *Trends Chem.*, **2020**, *2* (2), 92–106 Plutschack, et al., Chem. Rev., **2017**, *117* (18), 11796–11893





Photochemistry Reaction Scalability

Scaling flow photochemistry reactions

- Flow photochemistry strategies can help in the scale up of photocatalytic reactions:
 - We can increase the throughput by running the reactions for longer.
 - We can increase the number of reactors.
 - We can increase the radiant flux into the reaction.

Flow reactors have been shown to deliver a 150-fold increase in absorbed photon flux c.f. batch reactors

- As light (the photon) is the reagent in photochemistry reaction, by increasing the number of photons reaching the reaction mixture, it is possible to increase the reaction rate:
 - We can increase the light intensity into the same photochemistry reactor using the same platform to scale-up our reactions.

Moles photons per hour = λ (nm) × power (W) × 3×10⁻⁵

Number of moles of photons (einsteins) per hour at a given wavelength $(\lambda)^4$



Flow Photochemistry?

The benefits of flow chemistry are widely known; these translate to flow photochemistry with further benefits:

- Improved irradiation of the reaction:
 - Flow reactors improve light transmission to reaction.
- Improved reaction scalability:
 - Increasing light power enables more photons, increasing reaction rates and increasing throughput.
- Improved reaction selectivity:
 - Single wavelength LED light sources enable selectivity of reaction pathways.

- Improved reproducibility:
 - Controlled reaction exposure times with no over irradiation.
- Improved mixing and heat exchange:
 - Flow offers reproducible mixing and precise temperature control.
- Ability to perform multiphase chemistry:
 - Small reactor design allows ability to use gases, crucial in photochemical application.
- Increased safety:
 - Light proof reactors and power interlocks prevent exposure of user to high intensity light.







Summary and outlook

- Continuous photochemistry has still not reached maturity, especially at a production scale.
- Recent developments have driven interest in the synthetic and process chemistry communities.
- Further development will be reliant on multiple discipline; material scientists for light source technologies and chemical engineering in the design of reactors along with chemists to exploit this technique.
- The outlook is that Flow photochemistry will become a key technique in the Green Chemists toolbox.





Asynt Flow Photochemistry

Asynt offers a range of flow photochemistry reactors that cover a wide range of applications, reactor type with efficient temperature control and a range monochromatic LED light sources.









Asynt has range of modern flow photochemistry reactors for most scales, & knowledgeable support.



Please come and visit the Asynt booth on stand **H102**.



www.asynt.com

