

Advancing Reproducible Photochemistry

Why Photoreactor Design Matters in Modern Photocatalysis

1. Introduction

Photochemistry has become an increasingly important enabling technology across synthetic chemistry, medicinal chemistry, process development and materials science. As the use of photocatalytic methods continues to expand, researchers require reactor platforms capable of delivering reproducible irradiation, controlled thermal management and scalable workflows.

A recent academic benchmarking study from Imperial College London and BASF investigated the influence of photoreactor design on photon delivery and photochemical reaction performance. The work compared multiple commercially available photochemical reactor architectures using quantitative actinometry and representative photocatalytic transformations.

The study reinforced several important principles relevant to modern photochemistry:

- Efficient photon delivery is critical for reaction performance
- Reactor geometry strongly influences reaction kinetics
- Direct irradiation can improve photon transfer efficiency
- Temperature control and mixing consistency are important for reproducibility

These findings are highly relevant for laboratories seeking reproducible and scalable photochemistry workflows.

The Asynt Illumin8 photochemical reactor aligns closely with these design principles through its direct irradiation architecture, controlled reaction geometry and flexible workflow compatibility.

2. The Growing Importance of Photochemistry

Photochemistry has undergone significant growth in recent years due to advances in LED technology and the increasing adoption of photocatalysis in synthetic chemistry.

Modern photochemical methods enable:

- Mild reaction conditions
- Access to previously challenging transformations
- Improved selectivity
- Enhanced sustainability
- Efficient catalytic activation

Registered Address: Hall Barn Road Industrial Estate, Isleham, Cambridgeshire CB7 5RJ, UK



Applications now span a broad range of chemistry disciplines and industrial sectors.:

- Pharmaceutical discovery
- Medicinal chemistry
- Agrochemical development
- Process chemistry
- Materials science
- Academic research

As photochemical methods continue to move from exploratory academic workflows into routine synthetic and process chemistry applications, the need for reproducible and standardised reactor platforms is becoming increasingly important. Small variations in irradiation geometry, photon delivery and thermal management can lead to significant differences in reaction performance between laboratories. The publication reinforces the importance of carefully engineered reactor systems capable of delivering consistent and reproducible photochemical conditions.

3. Why Reactor Design Matters

Unlike conventional thermal chemistry, photochemical performance is heavily dependent on photon flux and irradiation geometry.

Photocatalytic reactions are fundamentally dependent on efficient excited-state generation. Inefficient photon transfer can limit catalyst activation, reduce reaction kinetics and introduce variability into photochemical workflows. Reactor designs capable of maximising photon utilisation can therefore significantly influence reaction efficiency and reproducibility

Key reactor design considerations include:

- Light intensity and wavelength
- Uniformity of irradiation
- Distance between light source and reaction mixture
- Optical path length
- Temperature management
- Stirring and mixing efficiency
- Scalability and workflow flexibility

As mentioned above, even relatively small changes in reactor architecture can significantly influence:

- Reaction rate
- Product conversion
- Catalyst efficiency
- Experimental reproducibility
- Scale-up performance



4. Direct Irradiation and Photon Delivery

One of the key observations from the academic benchmarking study was the importance of efficient photon transfer in photocatalytic workflows. The study compared direct irradiation reactor geometries with commercially available reflective reactor designs.

The authors concluded that direct irradiation systems provided substantially improved photon delivery efficiency compared with reflective architectures. The publication also noted that the Asynt Illumin8 direct irradiation reactor compared favourably during actinometric evaluation.

5. Quantitative Actinometry Benchmarking

Photon delivery was evaluated using ferrioxalate actinometry, a widely used quantitative technique for measuring photon flux in photochemical systems.





Figure 1. Ferrioxalate actinometry reaction

Irradiation Time	Reflective Reactor Design	Asynt Illumin8
10 seconds	8.71×10^{-6} mol	2.10×10^{-5} mol
20 seconds	1.63×10^{-5} mol	3.92×10^{-5} mol
30 seconds	2.43×10^{-5} mol	5.55×10^{-5} mol

Table 1 – Comparison of Average Photons Delivered over Irradiation time.

Across all irradiation intervals tested, the direct irradiation architecture delivered substantially higher measured photon flux.

The results also demonstrated highly linear irradiation behaviour, supporting stable and reproducible light delivery.



6. Reaction Benchmarking Highlights

The study investigated several representative photocatalytic transformations to evaluate how reactor design influence's reaction kinetics.

Trans-Stilbene Photoisomerisation:




Figure 2. trans-Stilbene photoisomerisation benchmark reaction

The photoisomerisation study demonstrated that direct irradiation resulted in:

- Faster conversion
- Improved early-stage reaction kinetics
- Shorter reaction completion times

Phenacyl Bromide Dehalogenation:




Figure 3. Phenyl bromide dehalogenation

The dehalogenation study demonstrated that direct irradiation resulted in:

- Significantly enhanced reaction rates under improved photon-delivery conditions
- Faster product formation
- More efficient photocatalytic performance

Nickel-Catalysed Cross-Coupling Amination:




Figure 4. Nickel-catalysed cross-coupling reaction

This more complex photocatalytic system highlighted the importance of controlled irradiation and consistent photon flux showing:

- Consistent photon delivery
- Stable irradiation conditions
- Temperature management
- Controlled reactor geometry

7. Designing for Reproducibility

One of the central themes throughout the publication was the importance of reproducibility in modern photochemistry.

The study highlighted several reactor characteristics that influence experimental consistency:

- Controlled irradiation geometry
- Consistent vial positioning
- Stable temperature management
- Uniform stirring performance
- Reliable photon delivery

These factors are increasingly important for:

- Reaction optimisation
- Method development
- High-throughput screening
- Cross-laboratory reproducibility
- Process development workflows

8. Illumin8: Supporting Reproducible Modern Photochemistry

The findings from the academic benchmarking study closely align with the engineering principles behind the Asynt Illumin8 photochemical reactor. The system has been developed to support reproducible photochemistry through controlled irradiation geometry, direct photon delivery and stable experimental conditions.

As modern photocatalytic workflows become increasingly sensitive to reproducibility, reactor architecture is no longer simply a hardware consideration - it is an important experimental variable influencing reaction performance, consistency and method transferability.

9. Practical Benefits for Modern Photochemistry Laboratories

The principles highlighted throughout the publication are highly relevant to day-to-day photochemistry workflows where reproducibility, consistency and experimental transferability are increasingly important.

Carefully controlled irradiation conditions can support:

- Faster and more reliable reaction optimisation
- Reduced experimental variability
- Improved confidence during method development
- More reproducible photocatalytic screening
- Greater consistency between users and laboratories
- Improved transferability of experimental conditions

10. Key Takeaways

- ✓ Efficient photon delivery strongly influences photochemical performance
- ✓ Reactor geometry plays a critical role in reaction kinetics and reproducibility

Registered Address: Hall Barn Road Industrial Estate, Isleham, Cambridgeshire CB7 5RJ, UK



- ✓ Direct irradiation architectures can improve photon transfer efficiency
- ✓ Consistent temperature management and mixing are important for reliable workflows
- ✓ Modern photochemistry increasingly requires scalable and reproducible reactor platforms

11. Conclusion

As photochemistry becomes increasingly integrated into mainstream synthetic workflows, reactor reproducibility and photon-delivery efficiency are likely to become critical experimental parameters rather than secondary equipment considerations.

The academic benchmarking study reinforced the importance of efficient photon delivery, controlled irradiation geometry and reproducible experimental conditions in modern photocatalysis. The results demonstrate that reactor architecture can significantly influence photon transfer efficiency, reaction kinetics and experimental reproducibility.

Importantly, the publication also highlights a broader challenge facing modern photochemistry: the need for standardised and reproducible reactor platforms capable of delivering reliable experimental performance across laboratories and applications.

The engineering principles behind the Asynt Illumin8 photochemical reactor have been developed specifically around these requirements, supporting modern photochemistry workflows where consistency, reproducibility and efficient photon delivery are essential.

As photochemical methods continue to expand across synthetic chemistry, medicinal chemistry and process development, carefully engineered reactor systems will play an increasingly important role in enabling reliable and reproducible photochemical science.



12. Learn More About Illumin8

Discover how controlled irradiation and reproducible reactor design can support modern photochemistry workflows.

Visit: <https://www.asynt.com/product/illumin8-parallel-photoreactor/>

Or contact the Asynt team to discuss your photochemistry applications.

13. References

Bailey, R. M.; Martínez-Aguirre, M.; Schaefer, B.; Crimmin, M. R.; Miller, P. W.
"An accessible and efficient 3D printed modular 'M-Arc' photoreactor"
Reaction Chemistry & Engineering (2026).

Actinometry comparison data and supplementary information were also reviewed during preparation of this white paper.

