

PRESS RELEASE

June 2007



Flexible New Controlled Lab Reactor System

Cambridge UK: The new ReactoMate from Asynt is a competitively priced laboratory reaction system that is available in convenient standard configurations or can be fully customised to meet customers' specific requirements.

Announcing the new product, Asynt MD Martyn Fordham said: "We became aware of a real need for a CLR (Controlled Lab Reactor) system with the flexibility to meet users' specific needs, but at a sensible price. Using the 35 years combined experience in the Asynt technical team and our access to the best technologies in laboratory synthesis, we have put together a highly competent yet competitively priced package that will handle all reaction volumes from 10ml to 30litres under precisely controlled conditions"

ReactoMate displays all the key features desirable in a CLR – precise control of reaction variables, a powerful direct drive overhead stirrer,

excellent visibility (the outer vacuum jacket option prevents ice formation at sub-ambient temperatures), quick release fittings for convenient vessel change, and the ability to operate under vacuum or an inert atmosphere. Comprehensive data logging facilities are also available for reaction variables.

E-mail sales@asynt.com call 01638 781709 or visit www.asynt.com for more information.

Asynt Ltd, based near Cambridge, was formed in July 2003 to develop, supply and support new and novel products and services for the organic chemist. The company developed the DrySyn range of products (launched in March 2004) and has successfully introduced partner company PharmaCore's custom synthesis services into the European market. Asynt products are developed by a team of chemists with 35 years combined experience and a mission to respond to the real needs of their colleagues in industry and academia. Asynt is committed to investment in new products and chemistries and plans to launch new ideas on a regular basis.

E-mail sales@asynt.com, call +44 (0)1 638 781 709 or visit www.asynt.com for more information.

For product enquiries please contact:
Asynt
Kathrin Hilpert +44 (0)1638 781709
Email: Kathrin@asynt.com

For media enquiries please contact:
Phoenix MarCom Ltd
Paul Carter + 44 (0) 1223 873318
Email: paul.carter@phoenixmarcom.co.uk