


...SMALL, SMART, SIMPLE, SILENT




**New compact fully automatic vacuum pump  
for rotary and parallel evaporation.**

CHEMISTRY VACUUM PUMP  
MD 1 C VARIO<sup>+</sup> \*

*\* Exclusively from Asynt*

## HIGH END RESULTS WITH ALL MAKES OF ROTARY EVAPORATOR, EVEN THE MOST BASIC

- Simple “One touch” vacuum control throughout the whole process, even for solvent mixtures
- Suitable for high boiling point solvents
- Oil-free diaphragm pump with outstanding chemical resistance
- Ultra low vibration and whisper quiet operation
- **Now with 2 mbar ultimate vacuum**



TECHNICAL DATA		MD 1C VARIO <sup>+</sup>
Vacuum controller		CVC 3000
Number of heads / stages		4 / 3
Max. pumping speed	m <sup>3</sup> /h	1.7
Ultimate vacuum (abs.)	mbar	2
Ultim. vac. (abs.) with gas ballast	mbar	4
Inlet / Outlet connection		Hose nozzle DN 10
Degree of protection		mm IP 20
Dimensions (L x W x H)	mm	263 x 193 x 299
Weight	kg	6.0

ORDERING INFORMATION	MD 1C VARIO <sup>+</sup>
100-120/200-230 V ~ 50-60 Hz UK	21696662

### ACCESSORIES

Rubber vacuum tubing DN 10 mm (686002)  
Air admittance valve VBM-B (674217)

BOOK YOUR FREE EVALUATION OF THE NEW MD 1C VARIO<sup>+</sup> TODAY



29 Hall Barn Road Industrial Estate, Isleham, Ely, Cambridgeshire, CB7 5RJ, UK

Tel.: +44 (0) 1638 781 309

Fax.: +44 (0) 1638 781 306

Further information at [www.asynt.com](http://www.asynt.com)