

The DrySyn Parallel Synthesis Kit & Starter Kit

No more oil baths or clamps needed!

The DrySyn Parallel Synthesis Kit allows the scientist to interchangeably heat and stir multiple reactions in any standard size tube or vial using any standard hotplate.

Complete kit includes: MULTI-E base plate to fit any hotplate stirrer, heat resistant handle and 1 pack of three standard reaction vial inserts in one size.

Starter kit includes: MULTI-E Base with heat resistant handle, and your choice of any three standard reaction vial inserts in any size from those detailed below.

Remember to specify which size inserts you require when ordering!



Popular purchase options

Complete kit part number: ADS3-PS-(size)
Starter kit part number: ADS3-PSS-(sizes)

Hole size (Ø mm)	Holes per insert	Suited for:
11.80 Ø	9	2 mL HPLC sample vials
15.00 Ø	5	1 dram vials
16.20 Ø	4	CEM microwave tubes
17.40 Ø	4	Biotage 2-5 mL microwave tubes & 2 dram vials
20.20 Ø	4	Robot tubes
20.40 Ø	4	24 mm boiling tubes
25.75 Ø	4	Wheaton reactor vials
27.65 Ø	4	20 mL Scintillation vials
28.20 Ø	4	Biotage 10-20 mL microwave tubes

Also Available

Part Number: ADS3-R-(size)

Pack of 3 x standard reaction vial inserts in **one size** from those detailed here.

Part Number: ADS3-R-A-(sizes)

Pack of 3 x standard reaction vial inserts in **any size** from those detailed here.

Part Number: ADS3-RT-17,20

Pack of 3 x tapered **17.20 mm** inserts for Biotage 0.5-2 mL 17.00 mm microwave vials

(Not currently available as part of the Parallel Synthesis kits)



Also available custom manufactured to your required size!