



Introducing... The DrySyn Parallel Synthesis Kit

No more oil baths or clamps needed!

The DrySyn Parallel Synthesis Kit allows the scientist to interchangeably heat and stir multiple reactions in any standard size tube or vial using any standard hotplate.

Kit Includes: Base plate to fit any hotplate stirrer, heat resistant handle and 1 pack of three standard reaction vial inserts in one size.

Don't forget to specify which size inserts you require when ordering!



ADS3-PS -(size)

Options:

ADS3-PS-11.80 (9 position)

Will fit 2ml HPLC sample vials

ADS3-PS-15.00 (5 position)

Will fit 1 dram vials

ADS3-PS-16.20 (4 position)

Will fit CEM microwave tubes

ADS3-PS-17.40 (4 position)

Will fit Biotage 2-5ml microwave tubes

ADS3-PS-20.20 (4 position)

Will fit robot tubes

ADS3-PS-24.40 (4 position)

Will fit 24mm boiling tubes

ADS3-PS-25.75 (4 position)

Will fit Wheaton reactor vials

ADS3-PS-27.40 (4 position)

Will fit 20ml Scintillation vials

ADS3-PS-28.20 (3 position)

Will fit Biotage 10-20ml microwave tubes

Now available: ADS3-RT-17.20

**Tapered 17.20mm inserts for Biotage
0.5-2ml 17.00mm microwave vials**

(Not currently available as part of this kit)



Effective. Sustainable. Safe. Brilliant.



enquiries@asynt.com
Tel: +44 1638 781709
www.asynt.com