

Introducing... The DrySyn Parallel Synthesis Kit

No more oil baths or clamps needed!

The DrySyn Parallel Synthesis Kit allows the scientist to interchangeably heat and stir multiple reactions in any standard size tube or vial using any standard hotplate.

Kit Includes: MULTI-E Base plate to fit any hotplate stirrer, heat resistant handle and 1 pack of three standard reaction vial inserts in one size.

Don't forget to specify which size inserts you require when ordering!



ADS3-PS -(size)

Options:

ADS3-PS-11.80 (9 position)

Will fit 2ml HPLC sample vials

ADS3-PS-15.00 (5 position)

Will fit 1 dram vials

ADS3-PS-16.20 (4 position)

Will fit CEM microwave tubes

ADS3-PS-17.40 (4 position)

Will fit Biotage 2-5ml microwave tubes

ADS3-PS-20.20 (4 position)

Will fit robot tubes

ADS3-PS-24.40 (4 position)

Will fit 24mm boiling tubes

ADS3-PS-25.75 (4 position)

Will fit Wheaton reactor vials

ADS3-PS-27.65 (4 position)

Will fit 20ml Scintillation vials

ADS3-PS-28.20 (4 position)

Will fit Biotage 10-20ml microwave tubes

Now available: ADS3-RT-17.20

**Tapered 17.20mm inserts for Biotage
0.5-2ml 17.00mm microwave vials**

(Not currently available as part of this kit)



Effective. Sustainable. Safe. Brilliant.

enquiries@asynt.com

Tel: +44 1638 781709

www.asynt.com

Introducing... The DrySyn Parallel Synthesis Starter Kit

No more oil baths or clamps

The DrySyn Parallel Synthesis Starter Kit allows the scientist to interchangeably heat and stir up to 27 parallel reactions using any standard hotplate. Choose any three inserts from our standard range to use with the DrySyn MULTI-E Base.

Kit Includes: MULTI-E Base with heat resistant handle, and your choice of any three standard reaction vial inserts in any size from those detailed below.

Don't forget to specify which size inserts you require when ordering!



ADS3-PSS -(size/size/size)

Options:

- | | |
|---|---|
| 11.80 (9 position) Will fit 2ml HPLC sample vials | 24.40 (4 position) Will fit 24mm boiling tubes |
| 15.00 (5 position) Will fit 1 dram vials | 25.75 (4 position) Will fit Wheaton reactor vials |
| 16.20 (4 position) Will fit CEM microwave tubes | 27.65 (4 position) Will fit 20ml Scintillation vials |
| 17.40 (4 position) Will fit Biotage 2-5ml microwave tubes | 28.20 (4 position) Will fit Biotage 10-20ml microwave tubes |
| 20.20 (4 position) Will fit robot tubes | Tapered 17.20mm inserts for Biotage 0.5-2ml 17.00mm microwave vials are not currently available as part of this kit. |



Effective. Sustainable. Safe. Brilliant.

enquiries@asynt.com
Tel: +44 1638 781709
www.asynt.com