

## Guide to choosing the correct magnetic stirring flea for your flask

### Optimise your stirring

The oval magnetic stirring bars are generally considered the most effective for stirring via a magnetic hotplate stirrer when using three flasks at a time via the DrySyn MULTI heating block system.



By choosing the correct size for your flask you can make a considerable difference to the quality and uniformity throughout your reactions.



### Popular purchase options

*Other sizes & shapes also available*

- ADS3-SB-15** Stirrer bar for 10ml flasks (pack of 3)
- ADS3-SB-25** Stirrer bar for 25ml and 50ml flasks (pack of 3)
- ADS3-SB-30** Stirrer bar for 100ml to 500ml flasks (pack of 3)



Contact us to discuss your requirements today