

CB Compact Scale

ADAM



Speed, Performance, Value - The Right Balance

- Easy to use keypad
- Stainless steel pan
- Portable with battery power
- g, oz and lb



www.adamequipment.com

Technical Specifications

Model	CB 501	CB 1001	CB 3000
Maximum capacity	500g	1000g	3000g
Readability	0.1g	0.1g	1g
Pan size	130mm & Oslash;		

General Specifications

Overall dimensions (wxdxh)	160 x 205 x 45mm
Net weight	0.385kg

Key Accessories

- In-use Cover
- Calibration Certificate
- M1 500g Calibration Weight for CB 501
- M1 1kg Calibration Weight for CB 1001
- 6vDC 200mA Adaptor



Adam Equipment - offices worldwide

www.adamequipment.com