

Racks

Tube storage made easier



Laboratory sample storage is made easier with Camlab's diverse range of racks suitable for microtubes, PCR tubes, test tubes and cryovials. The entire range has been developed to improve ways of working within a laboratory environment.

All our racks are manufactured from either polypropylene or ABS, enabling the racks to be stored in freezers with temperatures down to -80°C .

As a UK manufacturer of plastics, Camlab can offer bespoke products to help find the right solution for you.

All racks within the range are used for sample handling and preparation techniques on the benchtop.






For further details, samples and pricing please contact us using the information overleaf.



Features

- Racks available in a variety of colours for ease of organisation of samples.
- Fits into most freezer storage systems.
- Racks made from polypropylene, suitable for freezer storage from -80°C to 121°C .
- Racks made from ABS have high impact strength, good chemical resistance and mechanical toughness.
- Suitable temperature range from 40°C to 70°C .
- Clear lids available on 96 well reversible racks for sample enclosure.
- Stackable.
- Autoclavable.

Ordering Information

	Description	Part Number	Colour	Pack Size
	<p>48 and 96 Well Reversible Racks One side takes 0.1-0.5ml tubes flip over and the other side takes 1.5/2.0ml tubes. Designed with slightly larger wells to work easily with screw cap vials and ease of labelling. 96 place rack has a clear lid option available.</p> <p>Dimensions: 48 well rack (HxWxD) 21x68x198mm 96 well rack (HxWxD) 21x129x198mm</p>	ASY/7110-48B	Blue	Each
		ASY/7110-48G	Green	Each
		ASY/7110-48N	Natural	Each
		ASY/7110-48R	Red	Each
		ASY/7110-48Y	Yellow	Each
		ASY/7110-48BK	Black	Each
		ASY/7110-48MX	Mixed	Pk of 5
		ASY/7110-96B	Blue	Each
		ASY/7110-96G	Green	Each
		ASY/7110-96N	Natural	Each
		ASY/7110-96R	Red	Each
		ASY/7110-96Y	Yellow	Each
		ASY/7110-96MX	Mixed	Pk of 5
	<p>3D Benchtop Racks manufactured from ABS Variable positioning for 0.5 to 1.5ml micro tubes. Two rows of holes that accommodate two tubes per angle section. Holds up to 48 microtubes. Wipe-clean whiteboard for easy labelling. Supplied with marker pens. Dimensions: (HxWxD) 40x146x204mm</p>	ASY/7450-C	Cream	Each
		ASY/7450-P	Pink	Each
		ASY/7450-B	Blue	Each
	<p>81 Place Maxicold Freezer Racks Wells have lower conical profile to ensure microtubes remain upright. Lightweight and easily transportable, will fit into polycarbonate racking systems. Also available with slip on clear lid for multiple stacking. Both rack and lid have frosted area for writing on.</p> <p>Dimensions: (HxWxD) 27x137x145mm</p>	ASY/76000-B	Blue	Each
		ASY/76000-G	Green	Each
		ASY/76000-MX	Mixed	Pk of 5
		ASY/76000-N	Natural	Each
		ASY/76000-P	Pink	Each
		ASY/76000-Y	Yellow	Each
		ASY/76001-N	Natural + lid	Each
		ASY/76001-P	Pink + lid	Each
		ASY/76001-B	Blue + lid	Each
		ASY/76001-G	Green + lid	Each
		ASY/76001-Y	Yellow + lid	Each
		ASY/76001-MX	Mixed + lid	Pk of 5
		  	<p>3 Tier Autoclavable Tube Racks Convenient carrying handles on either end of the rack and ideal for use in a water bath as they will sink to the bottom can be stored at temperatures down to -80° C 60 place rack each hole is 16mm diameter. Low density and has a good strength to weight ratio. Offers excellent chemical resistance with very low permeability to water, vapor, and gases. Meets the requirement of the U.S. Food and Drug Administration Made from virgin polypropylene which avoids the presence of contaminating polymers, maintaining moulding consistency and ensuring plastic integrity. 3 Tier Rack Dimensions: (HxWxD) 62x135x275mm</p> <p>2 Tier Tube Racks for Ergonomically designed for use on the bench, in a water bath, shaker, fridge, freezer or incubator. For use with microtubes through to 50ml Falcon tubes. Base is manufactured in strengthened polypropylene and top section in ABS to improve impact resistance and toughness. Base tier has rounded wells with drain holes. Carry handle on end of the rack. 7020 racks hold up to 70 tubes, 7022 racks hold up to 18 falcon tubes NOTE: Racks are NOT autoclavable. 2 Tier Rack Dimensions: (HxWxD) 39x135x275mm</p>	ASY/7031-PB
ASY/7031-PBASS10	Pale Blue			Box 10
ASY/7031-PBFLAT15	Pale blue			Box of 15
ASY/7031-PK	Pale Pink			Each
ASY/7031-PKASS10	Pale Pink			Box 10
ASY/7031-PKASS15	Pale Pink			Box of 15
ASY/7031-PG	Pale Green			Each
ASY/7031-PGASS10	Pale Green			Box 10
ASY/7031-PGFLAT15	Pale Green			Box of 15
ASY/7031-PY	Pale Yellow			Each
ASY/7031-PYASS10	Pale Yellow			Box 10
ASY/7031-PYFLAT15	Pale Yellow			Box of 15
ASY/7020-R	Red			Each
ASY/7020-B	Blue			Each
ASY/7020-G	Green			Each
ASY/7022-R	Red			Each
ASY/7022-B	Blue			Each
ASY/7022-G	Green			Each
	<p>Polypropylene Floating Rack Holds up to 18 microtubes 1.5/2.0ml. Suitable for use in a waterbath or incubator. Handle in the centre of the rack for easy pick up from bath. Dimensions: (HxWxD) 40x90x133mm</p>	RTP/7120-W	White	Pk of 5