





## PCR Bubble

The PCR bubble is designed to provide a clean air environment for the preparation of PCR reagents. The PCR bubble is supplied as standard with a hepa filtered air inlet and UV decontamination lamp.

### Dimensions:



### Weight:

38 kg

### Materials of Construction:

Base	10mm HDPE - White, Black, Blue, Green or Red
Top	6mm Cast Clear Acrylic
Fan Filtration Module	White Powder Coated Mild Steel

### Utilities Requirements:

Power:	Europe:	1ph 220-240Volt AC
	USA:	1ph 90-115Volt AC

### Operational Detail:

Face Velocity	0.4m/s
Alarm Type	Hot Wire Anemometer
Lamp	LED
UV Lamp	36 WATT short wave uv decontamination

### Performance Specifications:

Air Cleanliness	Class 100
Filtration	H14 Hepa



+44 1638 781 709



+44 1638 781 706

[www.Lab-Bubble.com](http://www.Lab-Bubble.com)

[www.Asynt.com](http://www.Asynt.com)