

Technology for Vacuum Systems

Fast and controlled evaporation at the touch of a button

NEW CONTROLLER WITH INTUITIVE OPERATION

acuubran

PC 3001

Chemistry Pumping unit PC 3001 VARIO

PC 3001 VARIO

More time saving than before through simplified operation

• Intuitive Operation

through clear text menues and easy to handle rotary jog wheel control

- **Over 30% shorter process times** through continuous optimal high evaporation rates and speed control
- Ideal for high boiling point solvents with ultimate vacuum down to 2 mbar (even with gas ballast 4 mbar)

New vacuum controller with

intuitive jog wheel operation

• Safe reproducibility for complex applications via user programmable process cycles



Automatic boiling point detection

genuine fully automatic process control without user programming

• Proven Reliability

even in tough operating conditions (successor to the market leading PC 2001 VARIO)

TECHNICAL DATA		PC 3001 VARIO
Diaphragm pump		Chemistry diaphragm pump MD 1C VARIO-SP
Vacuum controller		CVC 3000
Max. pumping speed (ISO 21360) 50/60 Hz	m³/h	1,7
Ultimate vacuum (total)	mbar	2
Ultimate vacuum with gas ballast	mbar	4
Inlet connection (IN)		Hose nozzle DN 6/10 Hose nozzle DN 10
Outlet connection (EX) Cooling water connection		2 x hose nozzle DN 6
Interface		RS232C
Dimensions (L x W x H)	mm	294 x 328 x 480
Mass ca.	ka	7.7

Items supplied: Pumping unit with chemistry diaphragm pump MD 1C VARIO-SP mounted on pump support, vacuum controller CVC 3000 with integrated, ceramic pressure transducer and air admittance valve, exhaust vapour condensor, inlet separator, on/off switsch, cable with plug and instruction for use.

ORDERING INFORMATION		
200-230 V ~ 50-60 Hz	plug CEE	69 67 00
200-230 V ~ 50-60 Hz	plug CH	69 67 01
200-230 V ~ 50-60 Hz	plug UK	69 67 02
100-120 V ~ 50-60 Hz	plug US	69 67 03
ACCESSORIES		
Cooling water valve VKW-B, VAC	UU•BUS	67 60 13
Cable extension VACUU·BUS, 2 r	n	61 25 52
Y-adapter VACUU·BUS		63 66 56



Technology for Vacuum Systems

With 150 employees and over 40 years of experience, VACUUBRAND manufactures the most comprehensive range of laboratory and instrumentation vacuum pumps, gauges and controllers for rough and fine vacuum.

The product range comprises rotary-vane pumps, oil-free diaphragm pumps, complete pumping units, flexible vacuum systems and solutions for local area networks. The range is completed by a wide choice of fittings and accessories as well as vacuum gauges and controllers for rough and fine vacuum.

Our technical literature is only intended to inform our customers. The validity of general empirical values and results obtained under test conditions for specific applications depends upon a number of factors beyond our control. It is, therefore, strictly the users' responsibility to verify carefully the validity or suitability to their specific requirements. No claims arising from the information provided in this catalogue will consequently be entertained. Technical data are subject to change without notice.

Pictures may depict accessories which are not supplied as standard under the catalogue number printed.

© VACUUBRAND GMBH + CO KG 09/2007 · 99 91 84