






## Glove Bubble

The Glove Bubble is designed to provide a sealed environment for the control of atmospheres when handling sensitive samples. The Glove Bubble is supplied as standard with Gauntlet style Hypalon Gloves.

### Dimensions:



### Weight:

15 kg

### Materials of Construction:

Base	10mm HDPE - White, Black, Blue, Green or Red
Top	6mm Cast Clear Acrylic
Gloves	Hypalon Gauntlet

The glove bubble can be supplied with various gas or utilities connections. Simply supply or specify on order.