



PRESS RELEASE

Manifold Ensures Safe Secure & Efficient Reaction Cooling

The new **DrySyn water manifold** from **Asynt** optimises use of laboratory fume cupboard coolant use by converting a single water tap into three condensing streams and into one drain.

This novel lab accessory is equally effective linked to a recirculating chiller or a Liebig condenser and avoids daisy chaining multiple condensers, keeping everything in your fume cupboard space tidy.

Quick release connectors on the DrySyn water manifold enable rapid reaction set up/dismantling and form a seal on both sides when removing condenser tubing, diverting the coolant flow to the other connected condensers thus avoiding leaks. The water manifold is designed to be clamped to standard fume cupboard scaffolding and will keep your condenser tubing safe, secure and tidy. It can be rotated through 90 degrees to allow vertical or horizontal alignment.

To watch a video demonstrating how easy the DrySyn water manifold is to use please visit www.asynt.com/product/drysynwatermanifold/. For further information please contact Asynt on +44-1638-781709 or sales@asynt.com.

Asynt Ltd

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ
T: +44 (0)1638 781709 F: +44(0)1638 781706 sales@asynt.com www.asynt.com

Registered office: Eldo House, Kempson Way, Bury St Edmunds, Suffolk, IP32 7AR Registration No: 5160407
VAT No: GB 838 5592 82



Asynt is a leading supplier of affordable products, consumables and services for chemists in industry and academia. With staff of trained chemists - Asynt is able to draw upon this in-depth applications knowledge to provide a high level of customer support for its DrySyn Heating Blocks, Controlled Laboratory Reactors, Synthesis Tools, Evaporators, Circulators, Temperature Control Systems, Vacuum Pumps and Laboratory Safety Equipment.

OCTOBER 2014

asyntpr36.doc

Illustrative image: (image available on request)



For more information please contact:

Media: Dr Bill Bradbury

+44-208-546-0869 / info@primetek-solutions.com

Asynt Ltd

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ
T: +44 (0)1638 781709 F: +44(0)1638 781706 sales@asynt.com www.asynt.com

Registered office: Eldo House, Kempson Way, Bury St Edmunds, Suffolk, IP32 7AR Registration No: 5160407
VAT No: GB 838 5592 82