

NMR Tubes and Accessories

5 mm NMR Tube Boro 5.1 (ASTM Type 1 Class B glass)

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (µm)	Camber (µm)	Length (inch)	Packed in lots of
Boroeco-5-7	4,950 ± 0,030	0,43	25	75	7	100
Boroeco-5-8					8	
Boro300-7	4,950 ± 0,020	0,43	7	25	7	5
Boro300-8					8	
Boro400-7	4,950 ± 0,015	0,43	7	13	7	5
Boro400-8					8	
Boro500-7	4,950 ± 0,015	0,43	5	13	7	5
Boro500-8					8	
Boro600-7	4,950 ± 0,010	0,43	4	4	7	5
Boro600-8					8	
Boro800-7	4,950 ± 0,005	0,43	3	3	7	5
Boro800-8					8	



5 mm NMR Tube Boro 3.3 (ASTM Type 1 Class A glass)

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (µm)	Camber (µm)	Length (inch)	Packed in lots of
D200-5-7	4,950 ± 0,030	0,38	25	25	7	5
D200-5-8					8	
D200-5-9					9	
D300-5-7	4,950 ± 0,020	0,38	7	25	7	5
D300-5-8					8	
D300-5-9					9	
D400-5-7	4,950 ± 0,015	0,38	7	13	7	5
D400-5-8					8	
D400-5-9					9	
D500-5-7	4,950 ± 0,015	0,38	5	13	7	5
D500-5-8					8	
D500-5-9					9	
D600-5-7	4,950 ± 0,010	0,38	5	7	7	5
D600-5-8					8	
D600-5-9					9	




10 mm NMR Tube Boro 5.1 (ASTM Type 1 Class B glass)

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
Boroeco-10-7	9,95 \pm 0,03	0,60	254	51	7	100
Boroeco-10-8					8	
Boro300-10-7	9,98 \pm 0,02	0,60	76	38	7	3
Boro300-10-8					8	
Boro500-10-7	9,98 \pm 0,01	0,60	76	38	7	3
Boro500-10-8					8	



5 mm NMR Tube Quartz Glass

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
Qtz300-5-7	4,950 \pm 0,030	0,38	7	25	7	5
Qtz500-5-7	4,950 \pm 0,015	0,38	5	13	7	5



NMR Tube with Screw Cap (open or closed top)

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
ds300-5-7	4,950 \pm 0,020	0,38	7	25	7	5
ds300-5-7-S						
ds500-5-7	4,950 \pm 0,015	0,38	5	13	7	5
ds500-5-7-S						
Screwcap	closed top					
Screwcap-S	septum top					



NMR Tube with constriction

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
Boroct300-5-7	4,950 \pm 0,020	0,43	7	25	8	5
Boroct400-5-7	4,950 \pm 0,015	0,43	7	13	8	5
Boroct500-5-7	4,950 \pm 0,015	0,43	5	13	8	5



NMR Tube with PTFE Valve

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
ptfe300-5-7	4,950 \pm 0,020	0,38	7	25	7	1
ptfe300-5-8					8	
ptfe500-5-7	4,950 \pm 0,015	0,38	5	13	7	1
ptfe500-5-8					8	



NMR Coaxial Tube (Bruker = b, Varian = v)

Item number	O.D. (mm)	Wall thickness (mm)	Concentricity (μm)	Camber (μm)	Length (inch)	Packed in lots of
coaxialb-5-7	3,950 \pm 0,020	0,38	7	25	8	5
coaxialv-5-7	3,950 \pm 0,020	0,38	7	25	8	5

PTFE Liner (for 5 mm tubes)

Item number

DTL-5-4-7 4 mm



NMR Caps 5 mm

Item number

caps-5-100 100 pieces / package

caps-5-100-color

-a: aqua	-as: assorted	-bk: black
-b: blue	-f: fuchsia	-g: green
-o: orange	-pk: pink	-p: purple
-r: red	-s: sky	-w: white
-y: yellow		



PTFE Caps 5 mm

Item number

caps-5-10-ptfe 10 pieces / package

caps-5-100-ptfe

100 pieces / package



NMR Tube Rack

Item number

nmrrack for 72 tubes (5 mm)



NMR Tube Washer

Item number

washer ground joint NS 29 / 32

Pipettebrush

Item number

pipe170 170 mm 80 pieces / package

pipe300 300 mm 50 pieces / package



NMR Brush

Item number

brush 1 piece



Micro Torch

Item number	
torch	85 mm x 143 mm x 46 mm



Pasteur Pipettes

Item number		
pasteur150	150 mm	250 pieces / package
pasteur230	230 mm	250 pieces / package
pasteur310	310 mm	250 pieces / package



Pipette Bulbs

Item number		
bulbgrey	grey (2 ml)	1 piece
bulbyellow	yellow	100 pieces



Disposable needles with Luer tip

Item number		
needle08x120	0,8 x 120 mm	100 pieces / package
needle09x70	0,9 x 70 mm	100 pieces / package
needle08x40	0,8 x 40 mm	100 pieces / package



NMR Spatula

Item number		
spatel	18/8 stainless steel	3 x 215 mm

2-Part Disposable Syringes with Luer tip

Item number		
syr1	1 ml	100 pieces / package
syr2	2 ml	100 pieces / package
syr5	5 ml	100 pieces / package
syr10	10 ml	100 pieces / package
syr20	20 ml	100 pieces / package
syr50	50 ml	30 pieces / package



Rubber Ballon

Item number	
ballon	1 piece

For UK Sales and pricing please contact: www.asynt.com

sales@asynt.com

Telephone 01638 781709

Asynt