



PRESS RELEASE

How Choice of Reactor Vessel Affects Reaction Performance

Asynt has produced a **new report**** that compares the thermal performance and advantages of **vacuum jacketed** versus **non-vacuum jacketed reaction vessels**.

Data is presented for vacuum jacketed and non-vacuum jacketed reaction vessels connected to a Julabo Presto A30 temperature control system. The report demonstrates that a vacuum jacketed vessel not only achieved faster rates of cooling but was also able to maintain a reaction at -20°C compared to just -15°C for the non-vacuum jacketed vessel. Beneficially, also using a vacuum jacketed vessel was shown to eliminate ice formation on the outside of the reactor glass (a common problem with non-vacuum jacketed reaction vessels), allowing continuous viewing of the reactor contents.

By reducing heat loss / gain a glass vacuum jacket enables faster rates of heating and uses less energy to cool the system, thereby enabling lower temperatures to be achieved than reactor vessels without a vacuum jacket. The report concludes that using a vacuum jacketed reaction vessel will allow you to optimise thermal performance thereby improving control of reactions. In addition, the report proposes that the small additional cost of a vacuum jacketed vessel is easily justified not only by the additional performance but also by the cost saving of being able to use a lower powered temperature control system.

Asynt Ltd

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ
T: +44 (0)1638 781709 F: +44(0)1638 781706 sales@asynt.com www.asynt.com

Registered office: Eldo House, Kempson Way, Bury St Edmunds, Suffolk, IP32 7AR Registration No: 5160407
VAT No: GB 838 5592 82



Vacuum jacketed reaction vessels are available as options for the Asynt ReactoMate Pilot (5-50 litres) and for the Asynt ReactoMate Datum (100ml - 5 litres) reactor systems.

For a copy of the report and further information on vacuum jacketed reactor systems please visit <http://www.asynt.com/product/vacuum-jacketed-reaction-vessels/>** or contact Asynt on +44-1638-781709 / enquiries@asynt.com.

Asynt is a leading supplier of affordable products, consumables and services for chemists in industry and academia. With staff of trained chemists - Asynt is able to draw upon this in-depth applications knowledge to provide a high level of customer support for its DrySyn Heating Blocks, Controlled Lab Reactors, Synthesis Tools, Evaporators, Circulators, Temperature Control Systems, Vacuum Pumps and Laboratory Safety Equipment.

SEPTEMBER 2016

asyntpr62.doc

Asynt Ltd

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ
T: +44 (0)1638 781709 F:+44(0)1638 781706 sales@asynt.com www.asynt.com

Registered office: Eldo House, Kempson Way, Bury St Edmunds, Suffolk, IP32 7AR Registration No: 5160407
VAT No: GB 838 5592 82

Illustrative image: (image available on request)



For more information please contact:

Media: Dr Bill Bradbury

+44-208-546-0869 / info@primetek-solutions.com

Asynt Ltd

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ
T: +44 (0)1638 781709 F:+44(0)1638 781706 sales@asynt.com www.asynt.com

Registered office: Eldo House, Kempson Way, Bury St Edmunds, Suffolk, IP32 7AR Registration No: 5160407
VAT No: GB 838 5592 82