

DrySyn Classic

Single Position Heating Blocks

25 mL to 1000 mL

The original member of the DrySyn family, the DrySyn Classic allows the user to heat any size standard round-bottom flask from 25 mL to 1000 mL.

The base unit of the DrySyn Classic accepts a standard 1000 mL round bottom flask and allows the inserts to be used - one at a time - to adapt it to all standard flask sizes from 500 mL to 25 mL.

All single-position base units include two heat-resistant handles for safe handling when hot and adjustable feet to allow complete compatibility with any hotplate stirrer.



Ask us about other single position heating options!



DrySyn Classic

Single Position Heating Blocks

25 mL to 1000 mL

- ✓ Removes risk of oil fires
- ✓ No risks associated with oil spillages
- ✓ No toxic fume emission
- ✓ Reduction in electrical energy consumption
- ✓ No costly maintenance or disposal of oil



Popular Purchase Options*

ADS1: DrySyn Classic Kit

Includes 1000 mL base and 4 x adapter inserts for 50, 100, 250 & 500 mL round bottom flasks.


ADS1-25: DrySyn Classic 25 mL Adapter Insert

ADS1-A: DrySyn Classic Starter Kit

Includes 1000 mL base and 1 x adapter insert (specify size).

DrySyn Classic Heat Up Times

100 mL Flasks & Inserts, 50 mL Water



Ask us about other single position heating options!