



RAVF1825 Advanced Laboratory Freezer

Part Number	RAVF1825
Capacity (Litres)	505
Style	Freestanding
Exterior Dimensions (HxWxD)	1702 x 770 x 725mm
Interior Dimensions (HxWxD)	1490 x 645 x 545mm
Weight	86kg
Operating Temperature	-20°C
Interior Construction	White Enamel Steel
Exterior Construction	White Enamel Steel
Fixed Shelves	6
Door Type	Solid
Door Hinging	Right
Reverse Hinged	Optional
Defrost	Manual
Access Port	Optional
Open Door Alarm	X
High Temperature Alarm	✓
Low Temperature Alarm	✓
Remote Alarm Contacts	Optional
Digital Temperature Display	✓
Door Lock	✓
Feet	4
Energy Consumption	4.1 kWh/24h
Eco Friendly Refrigerant	✓

Labcold Advanced freezers are constructed from powder coated mild steel with an enamel steel interior so are built to last. Ideal for any application which requires the storage of a large quantity of product within tight temperature tolerances. This freezer is even manual defrost, putting you completely in control of the temperature of your samples.

**Looking for
temperature
mapping or
calibration?
Call the experts on
01256 705 570.**

Made to order at our UK factory, this freezer is specifically designed to withstand the rigours of professional use, and is backed by a comprehensive one year parts and labour warranty.

Fitted with a digital temperature display, you can be sure with a quick glance that the freezer is operating at the correct temperature. For extra peace of mind this freezer also has an audio/visual minimum and maximum temperature alarm and can have remote alarm contacts fitted at time of manufacture. For extra security this freezer comes complete with a door lock and two keys.