

# RLPR0214 Sparkfree Refrigerator



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chamber. Anything that could cause a spark is located outside of the chamber.

**Warning: Do not store open containers of volatile substances in this fridge.**



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm. The controller also records the minimum and maximum temperature since the last reset.

Part Number	RLPR0214
Capacity (Litres)	49
Style	Bench top
Exterior Dimensions (HxWxD)	565 x 480 x 455mm
Interior Dimensions (HxWxD)	420 x 400 x 320mm
Weight	19 kg
Energy Consumption	0.3 kWh/24hr
Operating Temperature	+0°C to +10°C
Factory Set Point	+4°C
Interior Construction	Moulded ABS
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right
Door Lock	✓
Door Open Alarm	✓
Shelves	2
Shelf dimension (WxD)	415 x 240 mm
Spare shelves available	✓
Feet	4
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Min/Max Recording	✓
Off Cycle Auto Defrost	✓
GWP Friendly Refrigerant	✓
2 Year Warranty (UK only)	✓