

RLPR1514 Sparkfree Fridge



The **Labcold** Sparkfree range is engineered to store flammable materials safely with no sources of ignition inside the chambers. Anything that could cause a spark is located outside of the chamber.

Warning: Do not store open containers of volatile substances in this fridge.



For extra security all sparkfree models are fitted with a built in lock and digital temperature display. The display features high and low audio and visual temperature alarms as standard in addition to a door open alarm. The controller also records the minimum and maximum temperature since the last reset.

Part Number	RLPR1514
Capacity (Litres)	400
Style	Free standing
Exterior Dimensions (HxWxD)	1855 x 600 x 585mm
Interior Dimensions (HxWxD)	1620 x 506 x 485mm
Weight	69kg
Energy Consumption	1.34kWh/24hr
Operating Temperature	0°C to +10°C
Factory Set Point	+4°C
Interior Construction	Moulded ABS
Exterior Construction	Painted Steel
Door Style	Solid
Door Hinging	Right (Can be reversed at time of manufacture)
Door Lock	✓
Door Open Alarm	✓
Adjustable Shelves	6
Shelf dimensions (WxD)	5 at 500 x 410mm 1 at 500 x 230mm
Feet	2 (adjustable)
Rollers	2 at rear
Temperature Display	✓
High Temperature Alarm	✓
Low Temperature Alarm	✓
Min/Max Recording	✓
Access Port	Optional Extra
Automatic Defrost	✓
Fan Assisted Cooling	✓