


# YOU'RE ON AIR ✨

WATERLESS DISTILLATION



## A CUSTOM DISTILLATION HEAD ALLOWS US TO FULLY ELIMINATE RUNNING WATER FROM ALL OUR TEACHING LABORATORY CONDENSING

Developed in partnership with Asynt (manufacturers of the CondenSyn air condenser), and made in-house by our glassblower, this is the first product of its kind and will be commercially available soon.

Air distillation in the lab contributes 15,000 litres to the water savings from the CondenSyn – now over 60,000 litres every year for our undergraduate practicals.

**All part of the SmartCondensing initiative – saving over 250,000 litres of waste water every year**

**WANT TO KNOW MORE? Find us @SustainableUEA**