

## PRESS RELEASE

### **Secure Clamping of Parallel Chemistry Apparatus**

Designed for chemists undertaking parallel syntheses the **Asynt DrySyn 3-Position Clamp** enables secure clamping of up to **three round bottom flasks** or **condensers** to a standard boss head.

With three sizes of the 3-Position Clamp available, to suit DrySyn MULTI-E, MULTI-M and MULTI-S heating block systems, these devices are also perfect to use with the Asynt CondenSyn range of waterless air condensers.

Each DrySyn 3-Position Clamp incorporates a novel fast grip / release mechanism that minimises set-up time, whilst the unique overall design enables the easy simultaneous removal of all three reactions from the heat source should the need arise.

Durably constructed from aluminium and stainless steel, with a replaceable rubber covering - the compact and affordably priced DrySyn 3-Position Clamp is available in a range of sizes to securely hold condensers and round bottom flasks from 5 mL to 500 mL.

#### **Asynt Ltd**

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ  
T: +44 (0)1638 781709 F: +44(0)1638 781706 enquiries@asynt.com www.asynt.com



For further information on the DrySyn 3-Position Clamp please visit <https://www.asynt.com/product/drysyn-multi-clamp/> or contact Asynt on +44-1638-781709 / [enquiries@asynt.com](mailto:enquiries@asynt.com).

Asynt is a leading supplier of affordable products, consumables and services for chemists in industry and academia. With staff of trained chemists - Asynt can draw upon this in-depth applications knowledge to provide a high level of customer support for its DrySyn Heating Blocks, Controlled Lab Reactors, Synthesis Tools, Evaporators, Circulators, Temperature Control Systems, Vacuum Pumps and Lab Safety Equipment.

**MARCH 2020**

**asyntpr105.doc**

**Asynt Ltd**

Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ  
T: +44 (0)1638 781709 F: +44(0)1638 781706 [enquiries@asynt.com](mailto:enquiries@asynt.com) [www.asynt.com](http://www.asynt.com)

Registration No: 5160407

VAT No: GB 838 5592 82

**Illustrative images:**



For more information please contact:

Media: Dr Bill Bradbury +44-208-546-0869 / [info@primetek-solutions.com](mailto:info@primetek-solutions.com)

**Asynt Ltd**  
Unit 29 Hall Barn Road Industrial Estate Isleham Cambridgeshire United Kingdom CB7 5RJ  
T: +44 (0)1638 781709 F: +44(0)1638 781706 [enquiries@asynt.com](mailto:enquiries@asynt.com) [www.asynt.com](http://www.asynt.com)

Registration No: 5160407 VAT No: GB 838 5592 82