

Overview: Huber KISS Series

May 2020



HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Cooling capacity kW	Pump performance				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
<a href="#">KISS K15</a>	-20...200	±0,05	2	0,2 @ 0°C	14	0,25	10,5	0,17	15	290 x 152 / 200	230V 1~ 50/60Hz
<a href="#">KISS K20</a>	-30...200	±0,05	2	0,35 @ 0°C	14	0,25	10,5	0,17	20	290 x 329 / 150	230V 1~ 50/60Hz
<a href="#">KISS K25</a>	-30...200	±0,05	2	0,35 @ 0°C	14	0,25	10,5	0,17	25	290 x 329 / 200	230V 1~ 50/60Hz
<a href="#">KISS K6</a>	-25...200	±0,05	2	0,15 @ 0°C	14	0,25	10,5	0,17	4,5	140 x 120 / 150	230V 1~ 50/60Hz
<a href="#">KISS K6s</a>	-25...200	±0,05	2	0,21 @ 0°C	14	0,25	10,5	0,17	4,5	140 x 120 / 150	230V 1~ 50/60Hz

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: [enquiries@asynt.com](mailto:enquiries@asynt.com)

Web: [www.asynt.com](http://www.asynt.com)