

Overview: Julabo CF Series
May 2020



Model	CF30	CF40	CF31	CF41
Working temperature range (°C)	-30 ... +150	-40 ... +150	-30 ... +200	-40 ... +200
Temperature stability (°C)	±0.03	±0.03	±0.02	±0.02
Cooling capacity (kW)	0.32	0.47	0.32	0.47
Heating capacity (kW)	2	2	2	2
Pump capacity flow rate (l/min)	15	15	22 ... 26	22 ... 26
Pump capacity flow pressure (bar)	0.35	0.35	0.4 ... 0.7	0.4 ... 0.7
Datasheet	CF30	CF40	CF31	CF41

Please contact us for further information.

Asynt Ltd
Tel: +44 (0)1638 781709
Email: enquiries@asynt.com
Web: www.asynt.com