

Overview: Julabo Dyneo Series  
May 2020



Model	DYNEO DD-200F	DYNEO DD-201F	DYNEO DD-300F	DYNEO DD-600F	DYNEO DD-601F	DYNEO DD-900F	DYNEO DD-1000F	DYNEO DD-1001F
Working temperature range (°C)	-20 ... +200	-20 ... +200	-25 ... +200	-35 ... +200	-35 ... +200	-38 ... +200	-50 ... 200	-38 ... +100
Temperature stability (°C)	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01
Cooling capacity (kW)	0.2		0.3	0.6	0.6	0.9	1	1
Heating capacity (kW)	2		2	2	2	2	2	1.5
Pump capacity flow rate (l/min)	8 ... 27		8 ... 27	8 ... 27	8 ... 27	8 ... 27	8 ... 27	8 ... 27
Pump capacity flow pressure (bar)	0.1 ... 0.7		0.1 ... 0.7	0.1 ... 0.7	0.1 ... 0.7	0.1 ... 0.7	0.1 ... 0.7	0.1 ... 0.7
Datasheet	<a href="#">DD-200F</a>	<a href="#">DD-201F</a>	<a href="#">DD-300F</a>	<a href="#">DD-600F</a>	<a href="#">DD-601F</a>	<a href="#">DD-900F</a>	<a href="#">DD-1000F</a>	<a href="#">DD-1001F</a>

Please contact us for further information.

**Asynt Ltd**

Tel: +44 (0)1638 781709

Email: [enquiries@asynt.com](mailto:enquiries@asynt.com)

Web: [www.asynt.com](http://www.asynt.com)