

**Overview: Julabo F-Series**  
May 2020



Model	F250	F500	F1000
Working temperature range (°C)	-10 ... +40	0... +40	0... +40
Temperature stability (°C)	±0.5	±0.5	±0.5
Cooling capacity (kW)	0.25	0.5	1
Pump capacity flow rate (l/min)	15	24	23
Pump capacity flow pressure (bar)	0.35	0.5	1
Datasheet	<a href="#">F250</a>	<a href="#">F500</a>	<a href="#">F1000</a>

**Please contact us for further information.**

**Asynt Ltd**

Tel: +44 (0)1638 781709

Email: [enquiries@asynt.com](mailto:enquiries@asynt.com)

Web: [www.asynt.com](http://www.asynt.com)