

Overview: Julabo FC-Series
May 2020



Model	FC600	FC1200	FC1600	FC600S	FC1200S	FC1600S	FC1200T	FC1600T	FCW2500T
Working temperature range (°C)	-20 ... +80	-20 ... +80	-20 ... +80	-10 ... +80	-15 ... +80	-15 ... +80	-10 ... +80	-15 ... +80	-25 ... +80
Temperature stability (°C)	±0.2	±0.2	±0.2	±0.2	±0.2	±0.2	±0.2	±0.2	±0.2
Cooling capacity (kW)	0.6	1.3	1.65	0.5	1.2	1.55	1.1	1.45	2.5
Heating capacity (kW)	1.2	1.2	1.2	1.2	1.2	1.2	1.2	1.2	1.2
Pump capacity flow rate (l/min)	20	20	20	22	22	22	28	28	28
Pump capacity flow pressure (bar)	0.5	0.5	0.5	1.2	1.2	0.2	3.5	3.5	3.5
Datasheet	FC600	FC1200	FC1600	FC600S	FC1200S	FC1600S	FC1200T	FC1600T	FCW2500T

Please contact us for further information.

Asynt Ltd
Tel: +44 (0)1638 781709
Email: enquiries@asynt.com
Web: www.asynt.com