

Model	F25-HE	F32-HE	F34-HE	FP50-HE	FPW50-HE	FP35-HL	F25-HL	F32-HL	F33-HL	FP50-HL	FPW50-HL
Working temperature range (°C)	-28 ... +200	-35 ... +200	-30 ... +150	-50 ... +200	-50 ... +200	-35 ... +150	-28 ... +200	-35 ... +200	-30 ... +200	-50 ... +200	-50 ... +200
Temperature stability (°C)	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01	±0.01
Cooling capacity (kW)	0.26	0.45	0.45	0.9	0.9	0.45	0.26	0.45	0.5	0.9	0.9
Heating capacity (kW)	2	2	2	2	2	2	2	2	2	2	2
Pump capacity flow rate (l/min)	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26	22 ... 26
Pump capacity flow pressure (bar)	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7	0.4 ... 0.7
Datasheet	F25-HE	F32-HE	F34-HE	FP50-HE	FPW50-HE	FP35-HL	F25-HL	F32-HL	F33-HL	FP50-HL	FPW50-HL

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com