

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 017T-H	-10...100	±0,2	2	1,7 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 020T-H	-20...100	±0,2	2	2 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 025T-H	-10...100	±0,2	2	2,5 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 040T-H	-10...100	±0,2	2	4 @ 15°C	26	3,0	--	--	3,5	--	400V 3~ 50Hz *
Unichiller 045T-H	-20...100	±0,2	2	4,5 @ 15°C	26	3,0	--	--	3,5	--	400V 3~ 50Hz *
Unichiller 055T-H	-10...100	±0,2	2	5,5 @ 15°C	57	5,6	--	--	5	--	400V 3~ 50Hz *
Unichiller 060T-H	-20...100	±0,2	2	6 @ 15°C	80	5,6	--	--	5	--	400V 3~ 50Hz *
Unichiller 070T-H	-10...100	±0,2	2	7 @ 15°C	84	5,6	--	--	5	--	400V 3~ 50Hz *
Unichiller 100T-H	-20...100	±0,2	2	10 @ 15°C	96	5,6	--	--	8,36	--	400V 3~ 50Hz *
Unichiller 110T-H	-10...100	±0,2	2	11 @ 15°C	96	5,6	--	--	8,36	--	400V 3~ 50Hz *
Unichiller 130T-H	-10...100	±0,2	2	13 @ 15°C	90	5,65	--	--	14	--	400V 3~ 50Hz *
Unichiller 160T-H	-10...100	±0,2	2	16 @ 15°C	96	5,6	--	--	14	--	400V 3~ 50Hz *
Unichiller 180T-H	-20...100	±0,2	2	18 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 200T-H	-20...100	±0,2	2	20 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 210T-H	-20...100	±0,2	2	21 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 250T-H	-20...100	±0,2	2	25 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 260T-H	-20...100	±0,2	2	26 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 300T-H	-20...100	±0,2	2	30 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 350T-H	-20...100	±0,2	2	35 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 017T	-10...40	±0,5	--	1,7 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 020T	-20...40	±0,5	--	2 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 025T	-10...40	±0,5	--	2,5 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 040T	-10...40	±0,5	--	4 @ 15°C	26	3,0	--	--	3,5	--	400V 3~ 50Hz *
Unichiller 045T	-20...40	±0,5	--	4,5 @ 15°C	26	3,0	--	--	3,5	--	400V 3~ 50Hz *
Unichiller 055T	-10...40	±0,5	--	5,5 @ 15°C	100	5,6	--	--	5	--	400V 3~ 50Hz *
Unichiller 060T	-20...40	±0,5	--	6 @ 15°C	80	5,6	--	--	5	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 070T	-10...40	±0,5	--	7 @ 15°C	84	5,6	--	--	5	--	400V 3~ 50Hz *
Unichiller 100T	-20...40	±0,5	--	10 @ 15°C	96	5,6	--	--	8,36	--	400V 3~ 50Hz *
Unichiller 110T	-10...40	±0,5	--	11 @ 15°C	90	5,6	--	--	8,36	--	400V 3~ 50Hz *
Unichiller 130T	-10...40	±0,5	--	13 @ 15°C	90	5,6	--	--	14	--	400V 3~ 50Hz *
Unichiller 160T	-10...40	±0,5	--	16 @ 15°C	96	5,6	--	--	14	--	400V 3~ 50Hz *
Unichiller 180T	-20...40	±0,5	--	18 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 200T	-20...40	±0,5	--	20 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 210T	-20...40	±0,5	--	21 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 250T	-20...40	±0,5	--	25 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 260T	-20...40	±0,5	--	26 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 300T	-20...40	±0,5	--	30 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 350T	-20...40	±0,5	--	35 @ 15°C	210	4,7	--	--	14	--	400V 3~ 50Hz *
Unichiller 017Tw-H	-10...100	±0,2	2	1,7 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 020Tw-H	-20...100	±0,2	2	2 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 025Tw-H	-10...100	±0,2	2	2,5 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 030Tw-H	-20...100	±0,2	2	3 @ 15°C	26	3,0	--	--	2,5	--	400V 3~ 50Hz *
Unichiller 040Tw-H	-10...100	±0,2	2	4 @ 15°C	26	3,0	--	--	2,5	--	400V 3~ 50Hz *
Unichiller 055Tw-H	-10...100	±0,2	2	5,5 @ 15°C	57	5,6	--	--	5,9	--	400V 3~ 50Hz *
Unichiller 060Tw-H	-20...100	±0,2	2	6 @ 15°C	80	5,6	--	--	5,9	--	400V 3~ 50Hz *
Unichiller 070Tw-H	-10...100	±0,2	2	7 @ 15°C	84	5,6	--	--	5,9	--	400V 3~ 50Hz *
Unichiller 100Tw-H	-20...100	±0,2	2	10 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 110Tw-H	-10...100	±0,2	2	11 @ 15°C	90	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 130Tw-H	-10...100	±0,2	2	13 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 160Tw-H	-10...100	±0,2	2	16 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 180Tw-H	-20...100	±0,2	2	18 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 200Tw-H	-20...100	±0,2	2	20 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 210Tw-H	-20...100	±0,2	2	21 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 220Tw-H12 CO2	-20...100	±0,2	12	22 @ 15°C	200	4,6	--	--	13	--	400V 3~ 50Hz *
Unichiller 250Tw-H	-20...100	±0,2	2	25 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 260Tw-H	-20...100	±0,2	2	26 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 300Tw-H	-20...100	±0,2	2	30 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 350Tw-H	-20...100	±0,2	2	35 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 500Tw-H	-10...100	±0,2	2	50 @ 15°C	220	4,7	--	--	17	--	400V 3~ 50Hz *
Unichiller 017Tw	-10...40	±0,5	--	1,7 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 020Tw	-20...40	±0,5	--	2 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 025Tw	-10...40	±0,5	--	2,5 @ 15°C	25	3,0	--	--	2,5	--	230V 1~ 50Hz
Unichiller 030Tw	-20...40	±0,5	--	3 @ 15°C	26	3,0	--	--	2,5	--	400V 3~ 50Hz *
Unichiller 040Tw	-10...40	±0,5	--	4 @ 15°C	26	3,0	--	--	2,5	--	400V 3~ 50Hz *
Unichiller 055Tw	-10...40	±0,5	--	5,5 @ 15°C	57	5,6	--	--	5,9	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 060Tw	-20...40	±0,5	--	6 @ 15°C	80	5,6	--	--	5,9	--	400V 3~ 50Hz *
Unichiller 070Tw	-10...40	±0,5	--	7 @ 15°C	84	5,6	--	--	5,9	--	400V 3~ 50Hz *
Unichiller 100Tw	-20...40	±0,5	--	10 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 110Tw	-10...40	±0,5	--	11 @ 15°C	90	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 130Tw	-10...40	±0,5	--	13 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 160Tw	-10...40	±0,5	--	16 @ 15°C	96	5,6	--	--	6,5	--	400V 3~ 50Hz *
Unichiller 180Tw	-20...40	±0,5	--	18 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 200Tw	-20...40	±0,5	--	20 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 210Tw	-20...40	±0,5	--	21 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 250Tw	-20...40	±0,5	--	25 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 260Tw	-20...40	±0,5	--	26 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 300Tw	-20...40	±0,5	--	30 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *
Unichiller 350Tw	-20...40	±0,5	--	35 @ 15°C	210	4,7	--	--	15	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	Capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unichiller 500Tw	-10...40	±0,5	--	50 @ 15°C	220	4,7	--	--	17	--	400V 3~ 50Hz *

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com