

Overview: Huber Unistat Tango Series
July 2020

HUBER Model	Temp. range °C	Temp. stability °C	Heating capacity kW	capacity kW	Cooling				Fill vol. Litres	Bath opening/ bath depth WxD / BD mm	Mains voltage
					Pressure l/min	bar	Suction l/min	bar			
Unistat Tango w	-45...250	±0,01	3	0,7 @ 0°C	55	0,9	--	--	1,5	--	230V 1~ 50Hz
Unistat Tango wl	-45...250	±0,01	3	0,7 @ 0°C	55	0,9	--	--	1,5	--	230V 1~ 50Hz

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com