

Overview: Julabo FT Series
July 2020

Model	<u>FT200</u>	<u>FT400</u>	<u>FT402</u>	<u>FT900</u>	<u>FT902</u>	<u>FT903</u>
Working temperature range (°C)	-20 ... +30	-40 ... +30	-40 ... +30	-90 ... +30	-90 ... +30	-90 ... +30
Temperature stability (°C)			±0.5		±1	±1
Cooling capacity (kW)	0.25	0.45	0.45	0.3	0.3	0.3

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com