

Overview: Julabo MAGIO Series
September 2020

Model	MAGIO MS-600F	MAGIO MS-601F	MAGIO MS-900F	MAGIO MS-1000F	MAGIO MS-310F
Order No.	9 032 704	9 032 705	9 032 706	9 032 707	9 032 713
Working temperature range (°C)	-35 ... +200.0	-35 ... +200.0	-38 ... +200.0	-50 ... +200.0	-30.0 ... +200.0
Temperature stability (°C)	+/-0.01	+/-0.01	+/-0.01	+/-0.01	+/-0.01
Cooling capacity (kW)	0.6	0.6	0.9	1	0.33
Heating capacity (kW)	2	2	2	2	2
Pump capacity flow rate (l/min)	16 ... 31	16 ... 31	16 ... 31	16 ... 31	16 ... 31
Pump capacity flow pressure (bar)	0.24 ... 0.92	0.24 ... 0.92	0.24 ... 0.92	0.24 ... 0.92	0.24 ... 0.92
Datasheet	MS-600F	MS-601F	MS-900F	MS-1000F	MS-310F

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com