

Overview: Julabo MAGIO Series
September 2020

Model	MAGIO MS-BC4	MAGIO MX-BC6	MAGIO MX-BC12	MAGIO MX-BC26
Order No.	9 032 504	9 033 506	9 033 512	9 033 526
Working temperature range (°C)	+20.0 ... +300.0	+20.0 ... +300.0	+20.0 ... +300.0	+20.0 ... +300.0
Temperature stability (°C)	+/-0.01	+/-0.01	+/-0.01	+/-0.01
Heating capacity (kW)	2	3	3	3
Pump capacity flow rate (l/min)	16 ... 31	16 ... 31	16 ... 31	16 ... 31
Pump capacity flow pressure (bar)	0.24 ... 0.92	0.24 ... 0.92	0.24 ... 0.92	0.24 ... 0.92
Datasheet	MS-BC4	MX-BC6	MX-BC12	MX-BC26

Please contact us for further information.

Asynt Ltd

Tel: +44 (0)1638 781709

Email: enquiries@asynt.com

Web: www.asynt.com