

Frozen in Time Ltd.

Manufacturers of Freeze Drying Machines
and Vacuum Cold traps

Lablyo Freezedrier

UK Manufactured

Bench top freeze drier with new Touch Screen controls. Data logging and display of temperature and vacuum on real time graph. Data downloadable via USB. 20cm diameter x 35cm deep ice condenser chamber with 5kg ice capacity. 3kg ice collection in 24 hours. Condenser temperature -55°C . Mobile temperature probe for product sampling.

Part No's

Lablyo	FIT/LYO/49/1
Vacuum pump for Lablyo	FIT/LYO/59/0

Can be used with the following freeze drying accessories:

- Acrylic drying chamber
- 6 height shelf rack
- 6 heated shelf rack
- 4 shelf vial closing device
- 8 port flask manifold
- 48 port ampoule manifold
- Spin Freezer



Asynt

Tel: +44 (0)1638 781709
Email: enquiries@asynt.com
Web: www.asynt.com