

Analytical & Precision Balances

enquiries@asynt.com

ENTRY LEVEL

Luna Analytical and Precision Balances



The balance perfectly suited for everyday lab use at an affordable price, providing users with easy to operate and reliable features without compromising on professional excellence.

- External and Internal calibration models available
- Readability: 0.1mg to 0.1g
- Capacity: 80g to 15000g

KEY FEATURES

- Large LCD display with 24mm high digits against a black background for clear, instant readings in any lighting conditions
- RS-232 and USB interfaces as standard
- Fully assembled draft shield

STANDARD

Solis Semi-Micro, Analytical and Precision Balances



The mid-range balance with high performance quality. Featuring outstanding precision and a high-resolution graphic display, Solis provides users with a range of integrated solutions for a variety of lab applications.

- External and Internal calibration models available
- Readability: 0.01mg to 0.01g
- Capacity: 120g to 8200g

KEY FEATURES

- Enhanced precision and readability for high accuracy weighing and counting
- Simply formulation with 99 ingredients presets
- RS-232 interface for send time-stamped data direct to printer or PC within Good Laboratory Practices (GLP) Guidelines

ADVANCED

Equinox Semi-Micro, Analytical and Precision Balances



Intuitive, touchscreen operation and advanced features that deliver remarkable precision to labs. Equinox is the error-free solution to challenging weighing tasks.

- External and Internal calibration models available
- Readability: 0.01mg to 0.01g
- Capacity: 120g to 8200g

KEY FEATURES

- 127mm-high touchscreen display with 800x480-pixel resolution
- RS-232 interface for send time-stamped data direct to printer or PC within Good Laboratory Practices (GLP) Guidelines
- 9 integrated applications ideal for highly detailed testing requirements: *Checkweighing, Statistical Analysis, Density Determination, Parts Counting, Percentage Weighing, Dynamic Animal Weighing, Formulation/Recipe, Weighing and Textile*

MOISTURE ANALYZERS

PMB Moisture Analyzers



Effortlessly record moisture test results. The Adam Equipment PMB provides quick, accurate and clear test results either on the production floor or out in the field. Automatic test-setting speeds up operation and the circular halogen heater lamp allows for even heating and drying of the sample.

- Model capacities: 50g - 200g
- Resolution: 0.02g - 0.01g

KEY FEATURES

- Suitable for grain, cosmetics and food testing
- Data can be transferred to a memory stick, printed, or sent to a PC
- Internal memory for up to 99 test results