

## CondensSyn MAXI for reflux or for distillation

*Details of the different configurations available*



**CondensSyn MAXI  
Reflux Inner**

**CondensSyn MAXI  
Outer**

**CondensSyn MAXI  
Distillation Inner**



The image (left) shows the CondensSyn MAXI with the distillation set up, as it would be used in conjunction with a lab reactor vessel.

The outer component is identical to that in the reflux setup, but the inner component has a small drain arm at the bottom that has a GL 18 thread on it.



### Popular purchase options

- GB-C-MAXI-350-P-B34 is a 350mm length in reflux configuration.
- GB-C-MAXI-350-S-B34 is a 350mm length in distillation configuration.
- GB-C-MAXI-450-P-B34 is a 450mm length in reflux configuration.
- GB-C-MAXI-450-S-B34 is a 450mm length in distillation configuration.

*Contact us now for further details & pricing*