

## DrySyn Multi Position Heating Blocks

*For vials and flasks up to 500 mL*

The DrySyn MULTI range of heating & cooling blocks enable clean, safe synthesis for multiple reactions in parallel. Developed by chemists, DrySyn is an established technology chosen by the leading laboratories around the world.

DrySyn allows for the complete replacement of unsafe oil baths in the lab. So whether you use round bottom flasks, tubes or vials there are DrySyn options available for all.



- ✓ Removes risk of oil fires
- ✓ No risks associated with oil spillages
- ✓ No toxic fume emission
- ✓ Reduction in electrical energy consumption
- ✓ No costly maintenance or disposal of oil
- ✓ Compatible with any standard hotplate stirrer

### *Flexible solutions*

With a wide range of fully-interchangeable inserts available for tubes, vials and round-bottom flasks, the DrySyn MULTI allows for fast, energy efficient heating. All MULTI bases include a heat-resistant handle for safe lifting.

By removing the need for multiple hotplate stirrers, DrySyn MULTI provides a cost effective, space-efficient solution.

Note: DrySyn MULTI rely on a wide pattern magnetic field for good coupling when using multiple flasks or tubes



**1 Reaction**



**3 Reactions**



**27 Reactions**

## Four bases. Endless possibilities.



**Uno**

Up to 1 x 500 mL



**MULTI-E**

Up to 3 x 100 mL



**MULTI-M**

Up to 3 x 250 mL



**MULTI-S**

Up to 3 x 500 mL



### Popular Purchase Options

*Further kits available plus individual components. Ask for details!*

**DrySyn MULTI-E Complete Kit** Includes MULTI-E base, 3 x 25 mL inserts, 3 x 50 mL inserts & 3 x 100 mL inserts.

**DrySyn MULTI-E Starter Kit** Includes MULTI-E base, 1 x 25 mL insert, 1 x 50 mL insert & 1 x 100 mL insert.

**DrySyn MULTI-M Complete Kit** Includes MULTI-M base, 3 x 100 mL inserts & 3 x 250 mL inserts.

**DrySyn MULTI-M Super Kit** Includes MULTI-M base, 3 x 25 mL inserts, 3 x 50 mL inserts, 3 x 100 mL inserts & 3 x 250 mL

**DrySyn MULTI-S Complete Kit** Includes MULTI-S base, 3 x 250 mL inserts & 3 x 500 mL inserts.

**DrySyn Parallel Synthesis Kit** Includes MULTI-E base and 3 x vial inserts of one size.

**DrySyn Parallel Synthesis Starter Kit** Includes MULTI-E base and 3 x vial inserts of different sizes.

**DrySyn MULTI-E / MULTI-M / MULTI-S 3 position clamps**

