



Provides a comprehensive range of production capabilities.

Losses are mitigated by the tank's low evaporation rates.

Efficiently stores samples in varying volumes.



The transfer unit facilitates easy transfer.

The capacity of each system has been optimized taking into account the evaporation rate.

Advanced insulation materials provide high thermal efficiency.

# Asynt

The storage tanks are manufactured and tested in accordance with the 2014/68/EU Pressure Equipment Directive.



MODEL	CT020	CT035	C050
Capacity (L)	20	35	50
Neck Diameter (mm)	57	57	57
Evaporation Rate (L/Day)	0.20	0.22	0.40
Material	Aluminum	Aluminum	Aluminum
Static Hold Time (Days)	159	125	100
Outer Dimensions (cm)	40 x 56,5	45,7 x 66	45,7 x 81,8
Mass (kg)	14	17	20

